

Name _____

Period _____

Partner _____

Date _____

Flame Tests

Prelab Questions

- 1) What is a spectator ion?
- 2) If you have a positive ion in solution do you also have to have a negative ion with it?
- 3) Why must the wires be extremely clean? What might happen if you touch them with your fingers?
- 4) If you use solutions containing both positive and negative ions (or of KNO_3 , NaNO_3 , and KCl), how can you tell which ions contribute to the color you see?

Procedure:

- 1) Wear safety goggles and follow all lab safety rules.
- 2) All of the solutions for this experiment are in bottle with platinum wires (inoculating loops) attached. You can share these with the other lab groups.
- 3) If necessary, clean the wire by dipping it in a test tube that contains 10 drops of 1M HCl (aq) and then heating the wire in the burner flame.
- 4) Heat the platinum wire in the hottest part of the burner flame until it glows but shows no color above the wire.
- 5) When the platinum wire is clean, dip the wire in the test tube containing LiCl (aq) solution and hold it in the hottest part of the burner flame.
- 6) Record your observations in Table 1 on the Report Sheet. Sometimes the color appears for a very short time so it may be necessary to dip the wire in the LiCl (aq) solution and put it in the flame several times.
- 7) Repeat for KCl (aq) and NaCl (aq), CuCl_2 (aq), CaCl_2 , SrCl_2 (aq), and BaCl_2 (aq).
- 8) Test the unknown solutions labeled A-E. Test the solutions the same way you did above and record your observations in Table 2 on the Report Sheet.
- 9) Before leaving the laboratory, wash your hands thoroughly with soap and water.

Data Table 1

Solution	Observations
LiCl (aq)	
KCl (aq)	
NaCl (aq)	
CuCl ₂ (aq)	
BaCl ₂ (aq)	
SrCl ₂ (aq)	
CaCl ₂ (aq)	

Data Table 2

Unknown A	
Unknown B	
Unknown C	
Unknown D	
Unknown E	

Post Lab Questions

- 1) All the solutions you tested were 0.5 M. How would the experimental results differ if the solutions were 0.10 M or 1.0 M?
- 2) What would you observe if you tested a solution that contained both LiCl and KCl?
- 3) If solutions of LiNO₃, KNO₃, and Na₂SO₄ are used, what different results would you expect?